

ALGER EQUIPMENT CO., INC.

320 Flightline

Lago Vista, TX 78645

Phone: 512-267-0383

FAX: 512-267-0452

TOLL FREE: 800-320-1043

[RETURN TO MAIN PAGE](#)

ALGER-LITE 220v. LOW VOLTAGE FOOTSWITCH INSTALLATION

1. Install ALGER-LITE cabinet on ceiling as per other instructions.

**** Caution ****

Do NOT lift ALGER-LITE by lite bar or put any undue force on lite bar as motor gears will be destroyed!!

2. Attach low voltage box to top of ALGER-LITE and hook up the seven wires from ALGER-LITE to low voltage box as per diagram.
3. Run the 4 ea. low voltage wires to the footswitch from the low voltage box as shown.
4. Connect 220V. input as shown to the two BLACK wires
5. The green wire from the footswitch is not connected using the 220V LOW VOLTAGE option.

