




Alger Equipment Co., Inc.
320 Flightline
Lago Vista, TX 78645


Ph 512.267.0383
800.320.1043
Fax 512.267.0452

Alger Light
with **LED**



Alger Light Positioning (Side View)

Ideally position the Alger Light with LED between the hip and knee if you work with the patient in the "reclined" position.



Alger Light Positioning (Overhead View)

Ideally center the Alger Light with LED over the center-line of chair. If there is obstruction in the center the light can be off-center approximately 12" to either side of center.

